

Měřené parametry



BIO TECHNOLOGY
bodiTECH

S: Serum P: Plasma WB: Whole Blood

C-Tip™ ▼ 10µL C-Tip™ 10µL ▼ 30µL C-Tip™ 30µL ▼ 50µL C-Tip™ 50µL

Product class	Test item	Sample type	Working Ranges
Cardiac	Tn-I plus	WB/S/P	0.01-15.00 ng/mL
	CK-MB	WB/S/P	3-100 ng/mL
	D-Dimer	WB/P	50-10,000 ng/mL
	NT-proBNP	WB/S/P	10-30,000 pg/mL
	Myoglobin	WB/S/P	2-500 ng/mL
Cancer	PSA	WB/S/P	0.5-100.0 ng/mL
	AFP	WB/S/P	5-350 ng/mL
	CEA	S/P	1-500 ng/mL
Diabetes	HbA1c ▼ 10µL	WB	(NGSP) 4-15 %
	Microalbumin	Urine	2-300 mg/L
Hormone	TSH ▼ 50µL	S/P	0.09-80.00 µIU/mL
	-	WB	0.4-80.0 µIU/mL
	T3	S/P	0.77-7.70 nmol/L
	T4	S/P	10.23-300.00 nmol/L
	FSH	S/P	1-100 mIU/mL
	Total βhCG	WB/S/P	5-50,000 mIU/mL
	LH	S/P	1-100 mIU/mL
	PRL	S/P	1-100 ng/mL
	Testosterone	S/P	1-10 ng/mL
	Cortisol	WB/S/P	80-800 nmol/L
Infection	CRP ▼ 10µL	WB/S/P	0.5-200.0 mg/mL
	PCT ▼ 50µL	WB/S/P	0.1-100.0 ng/mL
	PCT Plus ▼ 30µL	WB/S/P	0.02-50 ng/mL
	ROTA	Feces	1.3 x TCID ₅₀ /mL
	Rota / Adeno	Feces	qualitative
	NORO	Feces	GI 0.04 ng/mL GI 0.2 ng/mL
	HBsAg	WB/S/P	0-300 COI
	Anti-HBs	WB/S/P	0-500 mIU/mL
	Anti-HCV	WB/S/P	0-300 COI
	Dengue NS1 Ag	WB/S/P	1 x 10 ³ pfu/mL
	Dengue IgG/IgM	WB/S/P	qualitative
	Zika IgG/IgM	WB/S/P	qualitative
	HIV Ag/Ab	WB/S/P	P24: 2 IU/mL
Others	Ferritin	S/P	10-1,000 ng/mL