



Urine Sediment Analyzer Detergent

User Manual

【Product Name】Urine Sediment Analyzer Detergent

【Packaging Specification】300mL, 500mL

【Intended Use】Rinse the pipeline and flow cell

【Test Principle】Flow cell technique and high-speed imaging technique

【Main Ingredient】Sodium hypochlorite

【Storage Conditions & Shelf Life】

Stored at 2 °C~30°C in a dry place, sealed and protected from sunlight. Refer to the label for the shelf life. After opened, the shelf life is 30 days.

【For Instrument】Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Series Urinalysis System.

【Usage】See user manuals of Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Series Urinalysis System.

【Sample Requirement】Fresh urine sample without centrifugation.

【Performance Indicator】pH is 12.10 ± 0.5 at $(25 \pm 1) ^\circ\text{C}$

【Matters Needing Attention】

1. Avoid contact with skin and eyes;
2. Flush eyes with the flow of warm water for a few minutes if the product splashes into eyes (ensuring eyelids opened). Seek medical advice immediately if pain or swelling happens. Flush with soap water and water once skin contact occurs.

【Bibliography】

1. Ferris JA Comparison and standards of urine microscopy Lab Med 14:659-662, 1983
2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

【Symbol Explanation】

	Please read the User Manual		Stored at		Manufacturer
	Lot Number		Valid until		
	Conform to EU in vitro diagnostic medical equipment directives 98/79/EC				
	Authorised Representative		Catalogue number		

【Manual Approved and Modified Date】04/2019



DIRUI INDUSTRIAL CO., LTD.

95 Yunhe Street, New & High Tech.

Development Zone Changchun, Jilin 130012 P.R.China

Tel: +86 431 85100409 Fax: +86 431 85172581

E-mail: dirui@dirui.com.cn <http://www.dirui.com.cn>



EMERGO EUROPE Prinsessegracht 20
2514 AP The Hague The Netherlands