



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

Page 1 of 1

TO MEDISTA SPOL S.R.O.

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

INVOICE NO DATE OF ISSUE
1027734 06-AUG-25

CODE NO	DESCRIPTION	REF NO	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
814-0319	TRANSFER BAG 300ML 5DPC CE	-	250704001	1,000	04/07/25	JUL2028
814-0503	TRANSFER BAG 500ML CE	-	250704002	2,000	04/07/25	JUL2028
833-8441	BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF 16G (D)	-	250530001	72	30/05/25	MAY2028
833-8441	BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF 16G (D)	-	250703001	3,528	03/07/25	JUL2028

TEST ITEMS	JUDGEMENT	REMARKS
1.TOTAL APPEARANCE	PASSED	
2.FUNCTION TEST	PASSED	
3.PYROGENITY-USP	PASSED	
4.STERILITY-USP	PASSED	
5. SAFETY/TOXICITY-USP	PASSED	
6. CHEMICAL CONTENTS	PASSED	

TEST ITEMS 3-6 ARE NOT APPLICABLE TO BULK AND NON-STERILE PRODUCTS
TEST ITEM 6 IS APPLICABLE TO BLOOD BAG ONLY

P.O. NO : BE25/000131

CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 833-8441

LOT NUMBER : 250703001

DESCRIPTION : BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF
16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	100.8	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	99.9	PASSED
	* TOTAL CITRATE (%)	95 - 105	101.4	PASSED
	DEXTROSE (%)	95 - 105	99.8	PASSED
	TOTAL SODIUM (%)	95 - 105	99.7	PASSED

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.2	PASSED
	SODIUM (%)	95 - 105	99.7	PASSED
	CHLORIDE (%)	95 - 105	100.0	PASSED
	DEXTROSE (%)	95 - 105	99.2	PASSED
	MANNITOL (%)	95 - 105	101.5	PASSED
	ADENINE (%)	90 - 110	100.0	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 06-AUG-25

CY CHIA
QA MANAGER