



## Blood Collection Tube

Cell-Free DNA BCT® is a blood collection tube with a preservative stabilizing nucleated blood cells. This unique stabilization prevents the release of genomic DNA, allowing isolation of high-quality cell-free DNA. It also has been demonstrated to minimize the degradation of CTCs.

### Multiple Applications

Recover high-quality cell-free DNA suitable for a wide range of downstream applications commonly used in clinical research studies, drug discovery and diagnostic assay development.

### Consistent Results

Enhanced stability and recovery of cell-free plasma DNA minimizes variability associated with cell-free DNA sample preparation.

### Improved Logistics

Cell-Free DNA BCT reduces the need for immediate plasma preparation and CTC processing due to unique stabilization properties. Cell-free DNA is stable for up to 14 days at 6 °C to 37 °C, while circulating tumor cells are stable for up to 7 days, at room temperature, allowing convenient sample collection, transport and storage.

### Specifications

Blood Draw Volume	10.0 mL
cfDNA/gDNA Sample Stability	14 days at 6 °C to 37 °C
CTC Sample Stability	7 days at 15 °C to 30 °C

### Ordering Information

Cell-Free DNA BCT®	Catalog #
6-Tube Pack (10.0 mL) glass	218961
100-Tube Pack (10.0 mL) glass	218962
1000-Tube Pack (10.0 mL) glass	218992

Cell-Free DNA BCT® CE	Catalog#
6-Tube Pack (10.0 mL) glass	218996
100-Tube Pack (10.0 mL) glass	218997
1000-Tube Pack (10.0 mL) glass	230244

FOR RESEARCH USE ONLY. NOT FOR USE IN  
DIAGNOSTIC PROCEDURES IN THE U.S.

FOR EXPORT ONLY. NOT FOR SALE IN THE U.S.